

SHARPSGUARD®

Assembly



- 1 Match lid colour to label colour



- 2 'CLICK' lid onto container



- 3 Complete the label



If you have questions or need further information, please contact your local Area Manager or our Customer Service Team on +44 (0)1865 371841 or info.healthcareuk@mauserpackaging.com and visit our website www.daniels.co.uk.



REDEFINING SUSTAINABILITY

www.daniels.co.uk



info.healthcareUK@mauserpackaging.com

Unless otherwise stated, the information contained herein is at the specific request of the user and has been made available by Daniels Healthcare.

V02.2024.SC