

DANIELS® TRANSPORT CONTAINER RANGE ~ delivering packaging that conforms

The DANIELS® Transport Containers are a range of safe, easy to use and cost effective re-usable containers, designed to meet your transportation needs in the hospital, between healthcare sites and in the community.

The greatest feature of the Transport Container Range is that it helps you achieve compliance with the stringent performance test standards laid down by the United Nations and required by the Carriage of Dangerous Goods and Use of Transportable Pressure Equipment Regulations (2004.)

These regulations cover, for example, the transportation of clinical samples or specimens to a laboratory for diagnostic analysis.

Are you compliant with transport regulations?

- Packaging Instruction 650 or 621
- UN3373 or UN3291

DANIELS® Diagnostic Specimen Container 30

Nominal capacity 30 litres
Daniels product code DD854
NHS Supply Chain code KCP060


The 30 litre container is designed for larger collections from clinics, surgeries etc. to be carried on the public highway to the local hospital or laboratory for analysis. This is ideal for transport departments or carriers offering this type of collection service. The container has a double hinged lid which is secured in place by latches at the front corners. The container can be kept open to fill and then closed by pushing down on the lid ready for transportation.

DANIELS® Community Nursing Container (Large)

Daniels product code: DD855

The large Community Nursing Container is designed to be able to store large instruments that have been used at one location and then transported for cleaning at another. The container features a removable lid so that it can be reused.

The container is ideal for a midwife with lots of instruments to take to and from a home birth.



DANIELS® Diagnostic Specimen Container 7

Nominal capacity 7 litres
Daniels product code DD851
NHS Supply Chain code KCP043

This container is ideal as the outer packaging when transporting samples or specimens. The yellow label indicates the contents as Biological Substance - Category B.

Place the primary receptacle and the secondary packaging, (plastic bag), inside the container for a safe and effective means of transportation on the public highway. This works equally well for healthcare professionals in the community or hospital staff transporting samples to the laboratory. By utilising the same type of container in both hospital and community settings avoids any confusion of different practices.

The container has a fully removable lid, which is secured by four easily operated clips ready for transportation. Two removable dividers are provided to offer segregation of the contents. By using the in-car bracket provided the container can be securely stored in a car boot or footwell.

AVOID NON-COMPLIANCE!

Both the Diagnostic Specimen containers have been tested to UN drop test requirements and are compliant with Packaging Instruction 650, therefore, suitable for diagnostic specimens assigned to UN3373



DANIELS® Community Nursing Container (Small)

Daniels product code DD853
NHS Supply Chain code FSL262

This is an ideal container for community healthcare professionals who need a suitable container to hold small quantities of clinical waste, sharps and 'dirty' medical instruments, items used in the course of their everyday roles.

The container is supplied with two removable dividing inserts creating up to three separate compartments. These could hold both small sharps containers and/or small cytotoxic/cytostatic containers which will fit into the compartments at either end.

In this way the DANIELS® Community Nursing Container will assist with the segregation of waste as required by the 'Safe Management of Healthcare Waste' guidance issued by the Department of Health (2006).

The container has a fully removable lid, which is secured by four easily operated clips ready for transportation. By using the in-car bracket provided, the container can be securely stored in a car boot or footwell.

AVOID NON-COMPLIANCE!

The small Community Nursing container has been UN performance tested to meet all the requirements of Packaging Instruction 621 and is, therefore, suitable for clinical waste assigned to UN3291.

Daniels Healthcare is a market leader in the manufacture and supply of sharps and speciality containers.

If you would like to know more about this range, or any other Daniels Healthcare products, including the SHARPSGUARD® range, please call our Customer Support Team on 01865 371841



Technical Specification

	DANIELS® Diagnostic Specimen Container 7 Litre	DANIELS® Diagnostic Specimen Container 30 Litre	DANIELS® Community Nursing Container (Small)	DANIELS® Community Nursing Container (Large)
Order Code	DD851 <i>Also available from NHS Supply Chain Order Code KCP139</i>	DD854 <i>Also available from NHS Supply Chain Order Code KCP140</i>	DD853 <i>Also available from NHS Supply Chain Order Code FSL262</i>	DD855
Packaging	Carton of 10	Carton of 5	Carton of 10	Carton of 10
Nominal Capacity	7 Litres	30 Litres	7 Litres	22 Litres
Size (H x W x D, including lid)	240mm x 360mm x 130mm	355mm x 425mm x 305mm	240mm x 360mm x 130mm	345mm x 330mm x dia
Aperture Opening	146mm x 324mm	350mm x 266mm	146mm x 324mm	Removable lid
Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Quality Standards	UN Packaging Instruction 650 UN3373	UN Packaging Instruction 650 UN3373	UN Packaging Instruction 621 UN3291	UN Packaging Instruction 621 UN3291
Manufacturing Standards	ISO 9002	ISO 9002	ISO 9002	ISO 9002
Accessories	In-car bracket Daniels order code: DD507 (including with container)		In-car bracket Daniels order code: DD507 (including with container)	



Does your current packaging conform?



DH-QAD 195/25 Issue 7 23 July 2010

www.daniels.co.uk

DANIELS HEALTHCARE LTD
Unit 14, Station Field Ind. Est. Kidlington, OXON, OX5 1JD
Tel: 01865 371841



Safety, Service and Support - Every Day



Safety, Service and Support - Every Day